

# Kiyomi Healthcare India Pvt. Ltd.

## Kiyo-T<sup>®</sup>



Our Accreditations :



## Surface and Environment Disinfectant



# Kiyo-T<sup>®</sup> 125

### Composition

|  |        |
|--|--------|
| Alkyl Di Methyl Benzyl Ammonium chloride       | 2.37%  |
| Alkyl Di Methyl Ethyl Benzyl Ammonium Chloride | 2.37%  |
| Inert Ingredient                               | 95.26% |

Dilution : 15 ml. in 1 Lt. water for high risk area (for 1000 cft area)  
: 7.5 ml. in 1 Lt. water for low risk area (for 1000 cft area)

## Surface and Environment Disinfectant



# Kiyo-T<sup>®</sup> 256

### Composition

|   |         |
|---|---------|
| Didecyl Dimethyl Ammonium Chloride        | 8.70%   |
| N-Alkyl Dimethyl Benzyl Ammonium Chloride | 8.19%   |
| Insert ingredients                        | 83.11%  |
| Total                                     | 100.00% |

Dilution : 4 ml. in 1 Lt. water for high risk area (for 1000 cft area)  
: 2 ml. in 1 Lt. water for low risk area (for 1000 cft area)

## Instrument Disinfectant & Sterilization



# Kiyo-T<sup>®</sup> CDX

INSTRUMENT DISINFECTANT

### Composition

Glutaraldehyde -2.45% w/v with activator  
Inert ingredients

Contact Time : 30 - 45 Minute  
Long Lasting : Upto 28 Days Reuse life

## Instrument Disinfectant & Sterilization



# Kiyo-T<sup>®</sup> OPA

HIGH LEVEL INSTRUMENT DISINFECTANT

### Composition

Ortho- Phthaladehyde 0.55%  
Kiyo-T OPA is a high level disinfectant  
for surgical instrument  
and endoscope

Contact Time : 5 - 12 Minute  
Long Lasting : 14 Days Reuse life

## Antiseptic Microbicidal Soluton

# Kiyo-T<sup>®</sup> P.V



### Composition

Povidone Iodine - **Solution**  
5% & 10% w/v  
Povidone Iodine - **Scrub**  
7.5% w/v

Pack Size : 100, 500 ml & 2 Litre Bottle

## Antiseptic Microbicidal Ointment

# Kiyo-T<sup>®</sup> DON

Water Soluble Microbicidal Ointment



### Composition

Povidone Iodine IP : 5% w/w  
(Available Iodine : 0.5% w/w)  
Water Soluble Ointment Base : Q.S.

Pack Size : 250g Jar

## Surface and Environment Disinfectant



# Kiyo-T<sup>®</sup> Guard

### Composition

|  |        |
|--|--------|
| Alkyl Di Methyl Benzyl Ammonium chloride       | 2.37%  |
| Alkyl Di Methyl Ethyl Benzyl Ammonium Chloride | 2.37%  |
| Poly Hexamethylene Biguanide (PHMB)            | 2.00%  |
| Inert Ingredients                              | 93.26% |

Dilution : 15 ml. in 1 Lt. water for high risk area (for 1000 cft area)  
: 7.5 ml. in 1 Lt. water for low risk area (for 1000 cft area)

## Air & Environment Sterility Devices



# Kiyo-T<sup>®</sup> Fogg

STERILITY DEVICE

### Composition

Body & Tank : 304 Grade SS & Nylon Glass Filled  
Shock proof & Ultra Durable Body

Tank Capacity : 4 & 6 Ltr. (5 Ltr. & 7 Ltr. Tank)  
Blower Motor : 26000 RPM. Motor  
Range : 10 to 15 Feet  
Weight : 4.60 Kg. (Approx)

## Instrument Cleaning & Maintenance



# Kiyo-T<sup>®</sup> IR

RUST REMOVER

### Composition

Phosphoric Acid > 25% and Cleaner  
For cleaning of stains by immersion & ultrasonic cleaning of surgical instruments only made of stainless steel, hardened chrome steel or chrome-nickel steel.

Dilution : 10-100 ml in 1 ltr. Lukewarm water  
for 20 min. in immersion tray

## Instrument Cleaning & Maintenance



# Kiyo-T<sup>®</sup> Zyme

MULTI ENZYMATIC CLEANER

### Composition

Protease : Breaks peptide bond between amino acids & (Subtilisin) breaks down proteinaceous matter blood, protein rich body fluids  
Lipase : Breaks down body lipids and fatty substances, cream, ointments  
Amylase : Breaks down carbohydrates, polysaccharide, saliva & Mucus  
Cellulase : Breaks down cellulose, starchy matter, cell wall & cell membrane, faecal matter medicinal powder & liquids.

Dilution : 8 ml. in 1 Lt. water ( Immersion Tray )  
: 4 ml. in 1 Lt. water ( Ultrasonic )

## Antiseptic Microbicidal Soluton



# Kiyo-T<sup>®</sup> Lone

### Composition

|                            |              |
|----------------------------|--------------|
| Chlorhexidine Gluconate IP | : 0.3% (v/v) |
| & Cetrimide Solution       | : 0.6% (v/v) |
| Purified Water             | : q.s        |

Pack Size : 100 ml & 1 Litre Bottle

## Hemodialysis Disinfection



# Kiyo-T<sup>®</sup> Citro

DISINFECTION AND DECALCIFICATION OF HEMODIALYSIS MACHINE

### Composition


|                  |         |
|------------------|---------|
| Citric acid I.P. | 21% w/w |
| Lactic acid I.P. |         |
| Malic acid I.P.  |         |
| Purified water   | Q.S.    |

Consumption : Approx. 82/90/96ml per cycle

Pack Size : 5 Litre



C-43, DSIIDC Complex, Kalyanpuri, Delhi-110091

Customer Care No. : **93154 88324** 

Email : [kiyomihpl@gmail.com](mailto:kiyomihpl@gmail.com)

Web. : [www.kiyomihealthcare.in](http://www.kiyomihealthcare.in)

